

Note:

ClearNote / ClearNote⁺ / Light / VGA with Win8

ClearNote, ClearNote⁺, ClearNote Light and ClearNote VGA do work with Windows8.

For this please install the USB-Treibers STK1160Grabber with following steps:

- 1) From the Metro Start Screen, open Settings with press the keyboard Key "WIN" simultaneously with "C" (or move your mouse to the bottom-right-corner of the screen and wait for the pop-out bar to appear, then click the Gear icon).
- 2) Click 'More PC Settings'.
- 3) Click 'General'.
- 4) Scroll down, and click 'Restart now' under 'Advanced startup'.
- 5) Wait a bit.
- 6) Click 'Troubleshoot'.
- 7) Click 'Advanced Options'
- 8) Click 'Windows Startup Settings'
- 9) Click Restart.
Windows 8 start new and shows a menu
- 10) Press "F7"
- 11) Now you can install the STK1160